

Shamir Autograph Intelligence™

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.50	+0.75 to +3.50	Clear	11mm and up	15 S AI	S 250
Transitions®	1.5	-12.25 to +6.50	+0.75 to +3.50	Grey, Brown, Green & Style Colours*	11mm and up	TR S AI	S 250
Transitions® XTRActive™	1.5	-12.25 to +6.50	+0.75 to +3.50	Grey	11mm and up	X5 S AI	S 250
Polarized	1.498	-12.25 to +6.50	+0.75 to +3.50	Grey & Brown & Green	11mm and up	PZ S AI	S 250
DLC™ (Trivex)	1.53	-13.25 to +7.00	+0.75 to +3.50	Clear	11mm and up	TX S AI	S 250
DLC™ (Trivex) Transitions®	1.53	-13.25 to +7.00	+0.75 to +3.50	Grey & Brown	11mm and up	TTX S AI	S 250
DLC™ (Trivex) Blue Zero	1.53	-13.25 to +7.00	+0.75 to +3.50	Clear	11mm and up	BZX S AI	S 250
Polycarbonate	1.586	-14.75 to +8.75	+0.75 to +3.50	Clear	11mm and up	PC S AI	S 250
Polycarbonate Transitions®	1.586	-14.75 to +8.75	+0.75 to +3.50	Grey, Brown, Green & Style Colours*	11mm and up	PT S AI	S 250
Polycarbonate Transitions® XTRActive™	1.586	-14.75 to +8.75	+0.75 to +3.50	Grey	11mm and up	PX S AI	S 250
Polycarbonate Polarized	1.586	-14.75 to +8.75	+0.75 to +3.50	Grey & Brown & Green	11mm and up	PPZ S AI	S 250
Polycarbonate Blue Zero	1.586	-14.75 to +8.75	+0.75 to +3.50	Clear	11mm and up	BZP S AI	S 250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®

Shamir Autograph Intelligence™

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Plastic 1.6	1.6	-14.75 to +8.75	+0.75 to +3.50	Clear	11mm and up	16 S AI	S 250
Plastic 1.6 Transitions®	1.6	-14.75 to +8.75	+0.75 to +3.50	Grey & Brown & Green	11mm and up	T6 S AI	S 250
Plastic 1.6 Transitions® XTRActive™	1.6	-14.75 to +8.75	+0.75 to +3.50	Grey	11mm and up	X6 S AI	S 250
Plastic 1.6 Polarized	1.6	-14.75 to +8.75	+0.75 to +3.50	Grey & Brown	11mm and up	PZ6 S AI	S 250
Plastic 1.6 Blue Zero	1.6	-14.75 to +8.75	+0.75 to +3.50	Clear	11mm and up	6P S AI	S 250
Plastic 1.67	1.67	-16.75 to +10.00	+0.75 to +3.50	Clear	11mm and up	67 S AI	S 250
Plastic 1.67 Transitions®	1.67	-16.75 to +10.00	+0.75 to +3.50	Grey, Brown & Style Colours*	11mm and up	T67 S AI	S 250
Plastic 1.67 Transitions® XTRActive™	1.67	-16.75 to +10.00	+0.75 to +3.50	Grey	11mm and up	X67 S AI	S 250
Plastic 1.67 Polarized	1.67	-16.75 to +10.00	+0.75 to +3.50	Grey & Brown	11mm and up	PZ67 S AI	S 250
Plastic 1.67 Blue Zero	1.67	-16.75 to +10.00	+0.75 to +3.50	Clear	11mm and up	67P S AI	S 250
Plastic 1.74	1.74	-18.75 to +10.75	+0.75 to +3.50	Clear	11mm and up	74 S AI	S 250
Plastic 1.74 Transitions	1.74	-18.75 to +10.75	+0.75 to +3.50	Grey	11mm and up	74 S AI	S 250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®