

Shamir Autograph III®

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.00	0.75 to 4.00	Clear	11mm and up	15 S A ₃	S 250
Transitions®	1.5	-12.25 to +6.00	0.75 to 4.00	Grey, Brown, Green & Style Colours*	11mm and up	TR S A ₃	S 250
Transitions® XTRActive™	1.5	-12.25 to +6.00	0.75 to 4.00	Grey	11mm and up	X5 S A ₃	S 250
Polarized	1.498	-12.25 to +6.00	0.75 to 4.00	Grey & Brown & Green	11mm and up	PZ S A ₃	S 250
Drivewear®	1.5	-9.00 to +6.00	0.75 to 4.00	Drivewear®	11mm and up	DW S A ₃	S 250
DLC™ (Trivex)	1.53	-13.25 to +7.00	0.75 to 4.00	Clear	11mm and up	TX S A ₃	S 250
DLC™ (Trivex) Transitions®	1.53	-13.25 to +7.00	0.75 to 4.00	Grey & Brown	11mm and up	TTX S A ₃	S 250
DLC™ (Trivex) Blue Zero	1.53	-13.25 to +7.00	0.75 to 4.00	Clear	11mm and up	BZX S A ₃	S 250
Polycarbonate	1.586	-14.75 to +8.00	0.75 to 4.00	Clear	11mm and up	PC S A ₃	S 250
Polycarbonate Transitions®	1.586	-14.75 to +8.00	0.75 to 4.00	Grey, Brown, Green & Style Colours*	11mm and up	PT S A ₃	S 250
Polycarbonate Transitions® XTRActive™	1.586	-14.75 to +8.00	0.75 to 4.00	Grey	11mm and up	PX S A ₃	S 250
Polycarbonate Polarized	1.586	-14.75 to +8.00	0.75 to 4.00	Grey & Brown & Green	11mm and up	PPZ S A ₃	S 250
Polycarbonate Blue Zero	1.586	-14.75 to +8.00	0.75 to 4.00	Clear	11mm and up	BZP S A ₃	S 250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®

Shamir Autograph III®

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Plastic 1.6	1.6	-15.00 to +8.00	0.75 to 4.00	Clear	11mm and up	$\begin{matrix} 16 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.6 Transitions®	1.6	-15.00 to +8.00	0.75 to 4.00	Grey & Brown & Green	11mm and up	$\begin{matrix} T6 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.6 Transitions® XTRActive™	1.6	-15.00 to +8.00	0.75 to 4.00	Grey	11mm and up	$\begin{matrix} X6 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.6 Polarized	1.6	-15.00 to +8.00	0.75 to 4.00	Grey & Brown	11mm and up	$\begin{matrix} PZ6 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.6 Blue Zero	1.6	-15.00 to +8.00	0.75 to 4.00	Clear	11mm and up	$\begin{matrix} 6P \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.67	1.67	-16.75 to +9.00	0.75 to 4.00	Clear	11mm and up	$\begin{matrix} 67 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.67 Transitions®	1.67	-16.75 to +9.00	0.75 to 4.00	Grey, Brown & Style Colours*	11mm and up	$\begin{matrix} T67 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.67 Transitions® XTRActive™	1.67	-16.75 to +9.00	0.75 to 4.00	Grey	11mm and up	$\begin{matrix} X67 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.67 Polarized	1.67	-16.75 to +9.00	0.75 to 4.00	Grey & Brown	11mm and up	$\begin{matrix} PZ67 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.67 Blue Zero	1.67	-16.75 to +9.00	0.75 to 4.00	Clear	11mm and up	$\begin{matrix} 67P \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.74	1.74	-18.75 to +15.50	0.75 to 4.00	Clear	11mm and up	$\begin{matrix} 74 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$
Plastic 1.74 Transitions	1.74	-18.75 to +15.50	0.75 to 4.00	Grey	11mm and up	$\begin{matrix} 74 \\ \text{S} \\ A_3 \end{matrix}$	$\begin{matrix} \text{S} \\ 250 \end{matrix}$

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®