

Shamir Autograph II™

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.00	0.75 to 4.00	Clear	11mm and up	15 S 60	S 250
Transitions®	1.5	-12.25 to +6.00	0.75 to 4.00	Grey, Brown, Green & Style Colours*	11mm and up	TR S 60	S 250
Transitions® XTRActive™	1.5	-12.25 to +6.00	0.75 to 4.00	Grey	11mm and up	X5 S 60	S 250
Polarized	1.498	-12.25 to +6.00	0.75 to 4.00	Grey & Brown & Green	11mm and up	PZ S 60	S 250
DLC™ (Trivex)	1.53	-13.25 to +6.00	0.75 to 4.00	Clear	11mm and up	TX S 60	S 250
DLC™ (Trivex) Transitions®	1.53	-13.25 to +6.00	0.75 to 4.00	Grey & Brown	11mm and up	TTX S 60	S 250
DLC™ (Trivex) Blue Zero	1.53	-13.25 to +6.00	0.75 to 4.00	Clear	11mm and up	BZX S 60	S 250
Polycarbonate	1.586	-14.75 to +6.50	0.75 to 4.00	Clear	11mm and up	PC S 60	S 250
Polycarbonate Transitions®	1.586	-14.75 to +6.50	0.75 to 4.00	Grey, Brown, Green & Style Colours*	11mm and up	PT S 60	S 250
Polycarbonate Transitions® XTRActive™	1.586	-14.75 to +6.50	0.75 to 4.00	Grey	11mm and up	PX S 60	S 250
Polycarbonate Polarized	1.586	-14.75 to +6.50	0.75 to 4.00	Grey & Brown & Green	11mm and up	PPZ S 60	S 250
Polycarbonate Blue Zero	1.586	-14.75 to +6.50	0.75 to 4.00	Clear	11mm and up	BZP S 60	S 250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®

Shamir Autograph II™

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking	
Plastic 1.6	1.6	-15.00 to +6.00	0.75 to 4.00	Clear	11mm and up	¹⁶ S 60	S 250
Plastic 1.6 Transitions®	1.6	-15.00 to +6.00	0.75 to 4.00	Grey & Brown & Green	11mm and up	^{T6} S 60	S 250
Plastic 1.6 Transitions® XTRActive™	1.6	-15.00 to +6.00	0.75 to 4.00	Grey	11mm and up	^{X6} S 60	S 250
Plastic 1.6 Polarized	1.6	-15.00 to +6.00	0.75 to 4.00	Grey & Brown	11mm and up	^{PZ6} S 60	S 250
Plastic 1.6 Blue Zero	1.6	-15.00 to +6.00	0.75 to 4.00	Clear	11mm and up	^{6P} S 60	S 250
Plastic 1.67	1.67	-16.75 to +8.25	0.75 to 4.00	Clear	11mm and up	⁶⁷ S 60	S 250
Plastic 1.67 Transitions®	1.67	-16.75 to +8.25	0.75 to 4.00	Grey, Brown & Style Colours*	11mm and up	^{T67} S 60	S 250
Plastic 1.67 Transitions® XTRActive™	1.67	-16.75 to +8.25	0.75 to 4.00	Grey	11mm and up	^{X67} S 60	S 250
Plastic 1.67 Polarized	1.67	-16.75 to +8.25	0.75 to 4.00	Grey & Brown	11mm and up	^{PZ67} S 60	S 250
Plastic 1.67 Blue Zero	1.67	-16.75 to +8.25	0.75 to 4.00	Clear	11mm and up	^{67P} S 60	S 250
Plastic 1.74	1.74	-18.75 to +15.50	0.75 to 4.00	Clear	11mm and up	⁷⁴ S 60	S 250
Plastic 1.74 Transitions®	1.74	-18.75 to +15.50	0.75 to 4.00	Grey	11mm and up	⁷⁴ S 60	S 250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®