

Shamir Golf™

Product Availability and Markings

Material	Index	Power Range	Add	Color	Fitting Height	Invisible Marking
DLC™ (Trivex)	1.53	-5.00 to +4.00	1.00 to 3.50	Clear	19mm	TX  250
DLC™ (Trivex) Transitions®	1.53	-5.00 to +4.00	1.00 to 3.50	Grey & Brown	19mm	TTX  250
Polycarbonate	1.586	-5.50 to +5.25	1.00 to 3.50	Clear	19mm	PC  250
Polycarbonate Transitions®	1.586	-5.50 to +5.25	1.00 to 3.50	Grey, Brown, Green & Style Colours*	19mm	PT  250
Polycarbonate Transitions® XTRActive™	1.586	-5.50 to +5.25	1.00 to 3.50	Grey	19mm	PX  250
Polycarbonate Polarized	1.586	-5.50 to +5.25	1.00 to 3.50	Grey & Brown & Green	19mm	PPZ  250
Plastic 1.6	1.6	-6.50 to +5.25	1.00 to 3.50	Clear	19mm	16  250
Plastic 1.6 Transitions®	1.6	-6.50 to +5.25	1.00 to 3.50	Grey & Brown & Green	19mm	T6  250
Plastic 1.6 Transitions® XTRActive™	1.6	-6.50 to +5.25	1.00 to 3.50	Grey	19mm	X6  250
Plastic 1.6 Polarized	1.6	-6.50 to +5.25	1.00 to 3.50	Grey & Brown	19mm	PZ6  250

* Transitions Style Colours: Sapphire, Amethyst, Emerald Green, and Amber



SHAMIR

ReCreating Perfect Vision®