

NEW
from Shamir

Shamir FirstPAL™

NEW

Shamir FirstPAL™ –
Benefits for the patient,
Optometrist and Optician

- New unique optical solution for young presbyopic patients
- Special design makes it easier to enter the world of progressive lenses
- Maximum viewing comfort throughout the lens
- High patient satisfaction



SHAMIR



Shamir FirstPAL™

Finally, a solution for young presbyopes



Who?

Your patients in the young presbyopic 38-44 age range are beginning to feel a need for corrective lenses for close viewing and the majority of them may request a pair of reading glasses.

Why?

Young presbyopic patients are often not open to the idea of progressive lenses, mainly because they have heard that they may require a certain amount of adaptation.

The solution!

Now Shamir has created a lens specially designed to take the needs and concerns of this group into account. Shamir FirstPAL offers clear and wide far vision, a comfortable reading zone, and an easy transition through the intermediate zone. Easy on the eyes, Shamir FirstPAL provides the perfect solution for young presbyopic patients with a low addition power. Shamir FirstPAL gives them everything they need in one pair of glasses. It is affordable, and starting with a low addition it eases their introduction into the world of progressive lenses.

When ordering Shamir FirstPAL™ please provide the following information:

- Patient prescription (Progressive)
- Mono PD
- Fitting height
- Frame data: A, B and DBL

Shamir FirstPAL™, minimum fitting height 15mm

	RANGE	ADD POWER
1.5 Hard Resin Clear, Polarized, Transitions® VII & XTRActive®	-10.00 to +5.00	0.75 to 1.50
1.59 Polycarbonate Clear, Polarized, Transitions® VII & XTRActive®	-12.00 to +6.00	0.75 to 1.50
1.6 Plastic™ Clear, Blue Zero®, Polarized, Transitions® VII & XTRActive®	-12.00 to +6.00	0.75 to 1.50
1.67 Plastic™ Clear, Blue Zero®, Polarized, Transitions® VII & XTRActive®	-12.00 to +6.00	0.75 to 1.50



1828 Blue Heron Drive, Unit #37
London, ON N6H 0B7
www.shamir.com

